

**MITRA**

*Experience Peace of Mind*

# **M Stripper**

## **Cutter & Hand Sealer/Roller**



- **Stainless Steel Body**
- **Light Weight**
- **Spring Powered Handle**

# **M Stripper**

## **Cutter & Hand Sealer/Roller**

We are pleased to introduce MITRA 3-in-1 device in Hand Sealer, Cutter & Tube Stripper/Roller. As the name indicates, it is application for Dealing, Cutting and Squeezing the tubing filled with blood. Hand Sealer equipped with rollers for squeezing blood in tubing and seals tubing by aluminum ring.

- Completely stainless steel body, so no chance of rusting ever.
- Light weight
- Spring powered handle. To ensure uniform pressure while pressing to close and automatic recoiling of spring to release handle for opening.
- Screw-less rollers to avoid loosening of the rollers.
- Extra sharp cutting edges.
- Hot melt plastic coated handle covers to ensure no space or gap in between the handle arms and plastic covers.
- Ergometrically designed handle for better gripping.

### **TUBE STRIPPER/ROLLER**

- Keep the collection tube in between the rubber roller.
- Press the handle and roll it towards collection bag.

### **HAND SEALER**

- Fold the collection tube.
- Insert the Aluminum ring on it.
- Hold the folded tube with ring in up-down position and press between the grooves at the tip of the hand sealer.

FOR FURTHER INFORMATION PLEASE CONTACT:

**Mitra Industries (P) Ltd.**

A-181, Okhla Industrial Area, Phase-I, New Delhi - 110 020. INDIA.

**Tel.:** +91-11-26818971, 26818972, 26818973, 26810133, 41612545. **Fax:** +91-11-26810945, 26818970.

**e-mail:** enquiry@mipl.in **Internet:** www.mitraindustries.com