

MITRA*Experience Peace of Mind*

MSP

BLOOD COMPONENT BALANCE



MSP is a micro controller based system which measures the accurate volume & weight in ml & gm of the blood components according to the specific gravity of the component.

SPECIFICATIONS:

Load cell capacity	1200 grams
Weight sensor	Load cell
Weight display panel	Back light 4 line alphanumeric LCD display for better visibility
Indication Lamp	LCD display message LCD display in ml/gms
Calibration and Tarring	Auto
Accuracy	Weight calibration tolerance : ± 1 gm. Total accuracy of display : ± 1 gm.
Power source	220-240V AC 50Hz
weight	Approx 3 kg
Temperature characteristics	Operating : 0~30°C Storage : 0~40°C

Specification/Designs are subjected to change without any obligation on the part of the manufacturer

For Further Information Please Contact:

Mitra Industries (P) Ltd.

14/4, Delhi-Mathura Road, Near Metro Station, Mewala Maharajpur, Faridabad, Haryana- 121003, INDIA.

Ph.: 0129-2271068, 2271069. Toll Free No.: 1800 102 0307

E-mail: customersupport@mipl.in Website: www.mitraindustries.com