

MSP

BLOOD COMPONENT BALANCE



MSP is a micro controller based system which measures the accurate volume & weight in ml & gm of the blood components according to the specific gravity of the component.

SPECIFICATIONS:

Load cell capacity	1200 grams
Weight sensor	Load cell
Weight display panel	Back light 4 line alphanumeric LCD display for better visibility
Indication Lamp	LCD display message LCD display in ml/gms
Calibration and Tarring	Auto
Accuracy	Weight calibration tolerance : ± 1 gm. Total accuracy of display : ± 1 gm.
Power source	220-240V AC 50Hz
weight	Approx 3 kg
Temperature characteristics	Operating : 0–30°C Storage : 0–40°C

For Further Information Please Contact:

Mitra Industries Pvt. Ltd.

14/4, Delhi-Mathura Road, Near Metro Station, Mewala Maharajpur, Faridabad, Haryana - 121 003. INDIA

Phone: 0129-2271068, 2271069. Toll Free No.: 1800 1020 307

e-mail: customersupport@mipl.in Web: www.mitraindustries.com