



Give FREEDOM & PEACE to your PATIENT

- In other system during manufacturing , improper sealing between the two ports may lead to micro defects, but 4GYP has only one port.
- By virtue of their different contraction and expansion coefficient of the latex pillow (in green) and the surrounding PVC, gaps may arise between the two during the process of molding. An extra barrier (in red) in 4GYP prevents direct contact with the potential gaps around the latex and the PD solution as can occur in other connectology.
- The vacuum (in yellow) in the 4GYP port is an additional barrier for migration of infection into PD fluid.




Improve your
Quality of Life
With 4th Generation
CAPD SYSTEM



TRI-CHOICE

Multichamber Low GDP Bag for
Peritoneal Dialysis



CHOICE FOR LONG TERM PD...

- Preserves membrane and immune cell function
- Enhances Peritoneal Bio-compatibility
- No Infusion Pain
- Absolute patient comfort for those using Dialysis of more than one concentrations
- User friendly - Convenient to store and handle
- Ensures Long Life of Peritoneum
- Suitable for Long Term PD benefit



Design

- Portable & Compact
- Height adjustable stand for various patient position options



Advanced Technology

- Integrated heating system for cold climatic conditions
- Easy to use Touch Screen interface



Safety by Function

- Pre-attach tubing & drain bag (100% sterile & closed system)
- Unique emergency STOP, which re-starts from the stop sequence



Data Management

- USB flash drive compliant
(Transfer/visualize patient exchange data on desktop)



Automated Peritoneal Dialysis System

