

MSEAL

MITRA RF TUBE SEALER



M-seal BS-RF Tube Sealer is a bench-top instrument which gives uninterrupted sealing by switching over to battery support during sudden power failure. It gives more than 1200 seals. High sealing capacity makes it suitable for outdoor blood donation camps and in-house donations.

...extra features

- Automatic switch over to battery support during power failure.
- Light weight, portable has a compact design.
- With splash guard, provided to protect from splashing.
- Automatically triggers sealing with the help of sensor.
- LED indicator lamps for system ready, sealing & battery status.
- Short sealing time.
- Easy separation of tube system.
- No compressed air requirement.
- No warm up time.
- User friendly with ergonomical design.
- Carry box for the easy transportation.



...extra benefits

- m-seal easily creates hermetic snap-apart seals, user friendly.
- In-built battery charger.
- No RF cable required.
- Minimum charging time.
- Adapts to different tube sizes.
- Tube is easily pulled apart to separate seal ends.
- m-seal with sophisticated sensors which automatically adjust the sealing time as per the tube size and thickness.

Specifications

Input power supply	220V AC 50Hz 5A
Sealing suitability/applicability	Blood Bag tube
Sealing time	0.5-2 Sec.
Sealing capacity (battery)	Up to 1200*
RF output frequency	40.68 Mhz
Battery type	14.8 V - Li-ion
Inbuilt fuse	4A main supply AC
Mechanical size	W=3.4" H=8.2" L=12.9"
Weight	Approx = 4.5kg (with battery)



For Further Information Please Contact:

Mitra Industries Pvt. Ltd.

14/4, Delhi-Mathura Road, Near Metro Station, Mewala Maharajpur, Faridabad, Haryana - 121 003. INDIA

Phone: 0129-2271068, 2271069. Toll Free No.: 1800 1020 307

e-mail: customersupport@mipl.in Web: www.mitraindustries.com