




HONTO BILIARY STENT



Covered Biliary Stent



Uncovered Biliary Stent



Uncovered Biliary Stent for PTBD

Indicated for use in the palliative treatment of biliary strictures produced by malignant neoplasms.

Features:

HONTO Covered Biliary Stent

- Less invasive round edges.
- Excellent flexibility to fit naturally in bile duct.
- Optimum radial force.
- Excellent visibility in fluoroscopy.

HONTO Uncovered Biliary Stent

- Highly close cell construction.
- Excellent visibility under fluoroscopy.
- Interlock design.
- Excellent flexibility to fit naturally in tract.

HONTO Uncovered Biliary Stent

Product Name	Stent		Delivery Device	
	Diameter (mm)	Length (mm)	Diameter (Fr)	Usable Length(cm)
TBL10060180C0ZST	10	60	7.2	180
TBL10080180C0ZST	10	80	7.2	180
TBL10100180C0ZST	10	100	7.2	180

HONTO Covered Biliary Stent

Product Name	Stent		Delivery Device	
	Diameter (mm)	Length (mm)	Diameter (Fr)	Usable Length(cm)
TBL10060180C2XST	10	60	8.5	180
TBL10080180C2XST	10	80	8.5	180
TBL10100180C2XST	10	100	8.5	180

HONTO Uncovered Biliary Stent for PTBD

Product Name	Stent		Delivery Device	
	Diameter (mm)	Length (mm)	Diameter (Fr)	Usable Length(cm)
TBL10060060C0ZPT	10	60	7.2	60
TBL10080060C0ZPT	10	80	7.2	60
TBL10100060C0ZPT	10	100	7.2	60
TBL10120060C0ZPT	10	120	7.2	60

*Specification/Designs are subject to change without any obligation on the part of the manufacturer.



Manufactured & Marketed by:

Mitra Industries Private Limited

14/4, Delhi-Mathura Road, Faridabad, Haryana-121 003. INDIA

T: +91-129-4099300 | Toll Free: 1800 1020 307

E: customersupport@mipl.in | W: www.mitraindustries.com

CT-ST03, Rev.00
MIPL-MAR/2021