

TWIN BAG

Your first choice for long term Dialysis with CAPD



In CAPD, the dialysate flows into abdomen and stays there for a prescribed period of time (dwell time) usually 4-6 hours. Dextrose in the dialysate helps filter waste, chemicals and extra fluid in your blood from tiny blood vessels in the lining of your abdominal cavity.



Transfer Set



Transfer Set Cap



Titanium Adapter



Clamp

Advantages:

- Dialysis at home, or at work.
- Flexi hours 24 hrs. a day and 7 days a week.
- Exchange manually, 3-4 times a day, each exchange will take 25-30 minutes.
- Save travel time and money.
- User friendly – Convenient to store and handle.
- Ensure longer and happy life.
- Suitable for all ages.

Accessories required for CAPD

- Transfer Set
- Transfer Set Cap
- Titanium Adapter
- Clamps

Packing and Storage details

- CAPD bags should be stored between 4 to 30°C in a clean, dry, free from pests and rodents area.
- The bags should be protected from excessive heat.
- CAPD bags are available in volumes:- 1ltr, 2 ltr, 2.5 ltr, 5 ltr

Available dextrose concentration:- 1.5%, 2.5%, 4.25%.

Packing:- 1x10 = 10 bags/carton (1ltr)

1x6 = 06 bags/carton (2ltr)

1x5 = 05 bags/carton (2.5ltr)

1x3 = 03 bags/carton (5ltr)

These bags are packed in corrugated cartons as per above packing.

MITRA CAPD bags are available across India
Your first choice for excellent results

**CAPD is to be used on prescription of a registered medical practitioner only.*



Manufactured & Marketed by:

Mitra Industries Private Limited

14/4, Delhi-Mathura Road, Faridabad, Haryana-121 003. INDIA

T: +91-129-4099300 | Toll Free: 1800 1020 307

E: customersupport@mipl.in | W: www.mitraindustries.com