


ENDOPASS

Single use Colonoscope Distal Attachment Device




Device To Improve Colon Polyp Detection.
Worlds Only Graduated Mobility Ribbed Wing Design.

Features:

1. Colonoscopy is widely considered to be the best screening modality for colon cancer and adenomatous polyps by most gastroenterologists.
2. Significant number of Polyps are missed during colonoscopy screening and missed pool rate being reported to be 22%.
3. Some polyps maybe located on the proximal aspect of colon folds and maybe difficult to visualize during colonoscopy, especially during scope withdrawal.
4. This device is designed to improve the colon polyp detection by increasing mucosal surface area of visualization.

Specifications: Endopass (Colonoscope Distal Attachment)

| Model No. | Bush ID (mm) | Bush OD (mm) | Bush Length (mm) |
|------------------|--------------|--------------|------------------|
| T03CDA13M17ML27S | 13 | 17 | 27 |
| T03CDP11M15ML27S | 11 | 15 | 27 |



*Specification/Designs are subject to change without any obligation on the part of the manufacturer.



Manufactured & Marketed by:
Mitra Industries Private Limited
14/4, Delhi-Mathura Road, Faridabad, Haryana-121 003. INDIA
T: +91-129-4099300 | Toll Free: 1800 1020 307
E: customersupport@mipl.in | W: www.mitraindustries.com

CT-ES01, Rev.00
MIPL-MAR/2021